



prevention • monitoring • therapy • professiona

3in1 Respiratory Solution

Model: NK1000







Rossmax 3BREATH is a unique 3in1 compressor nebulizer that integrates with a nasal washer and a nasal aspirator. The nebulizer enables efficient medication delivered to the upper airways and lung periphery. It provides comprehensive, effective and convenient care for common respiratory conditions.

Piston Nebulizer

The Valve Adjustable Technology (VAT) bottle allows users to adjust different levels of Nebulization rate 0.1 (closed) / 0.3 (fully open) ml/min at consistent particle size.

- MMAD ≤2.5 µm Fine Particle Dose (FPD): 80-85%
- Patented Valve Adjustable nebulizer bottle (VAT)
- Very low noise with silencer and vibration absorber design
- Storage and Transportation Condition:
 -20°C~60°C(14°F~140°F), 10%~90% RH, 700~1060 hPa

Nasal Washer

- Flushes nasal cavities and rhinopharynx to clear secretions, mucus and harmful substances and alleviates allergic symptoms.
- Nasal washer nebulizes 10ml of solution in about 2 minutes.
- Silicone nasal adaptors provide comfortable and soft touch and reliable seal.
- Made of biocompatible and non-allergenic materials.
- Comes with 3 sizes of nasal adaptors and suits all ages.
- Fluid chamber capacity: Max 20ml

3 Nasal Aspirator

 Removes mucus to alleviate nasal congestion by clearing a blocked nose.
 It is safe and gentle solution to a blocked nose in babies.

Thanks to the Dual pump design, nebulization and suction pump, provides the maximum hygiene to avoid cross infection.

- Adjustable vacuum level.
- Ergonomic design to ensure a smooth operation.
- Max vacuum @ pump inlet: ≤ -450 mmHg
- Max air flow @ pump inlet: ≥ 5.0 lpm

Model	Qty per carton	Carton volume	QTY per 20'	QTY per 40'	QTY per 40'HQ
NK1000	4 pcs	0.042 cbm/ctn	2664 pcs	5300	6140 pcs