

Nebulizer Accessories Pack



N1 Patent Nos :
ZL 2012 2 0443597.3



Valve Adjustable Technology

Controllable nebulization rate at consistent particle size

The proprietary adjustable valve is able to deliver medications of different viscosity level according to every user's conditions without exchanging parts.

The Valve Adjustable Technology (VAT) bottle allows users to adjust different levels of Nebulization rate 0.08 (closed) / 0.5 (fully opened) ml/min at consistent particle size. Higher nebulization rate (fully opened) is for higher viscosity medications and higher breathing capacity user while lower nebulization rate with closed valve will be more appropriate for kids/infants with lower breathing capacity.

VAT is recommended for respiratory disorders such as Chronic Obstructive Pulmonary Disease (COPD) and asthma as well as intensive use in general.



- NB60/NH60/NI60: Fully open Valve > 0.2ml/min | Closed Valve > 0.08ml/min
- NB80: Fully open Valve > 0.3ml/min | Closed Valve > 0.1ml/min
- NA100: Fully open Valve > 0.4ml/min | Closed Valve > 0.15ml/min
- NB500: Fully open Valve > 0.5ml/min | Closed Valve > 0.15ml/min



Mouthpiece



Nebulizer kit (N1)



Filter
(5 pieces)



Air tube
(210cm)

Adult mask



Child mask

Model	Qty per carton	Carton volume
Neb Accessories Pack	50 pcs	0.082 cbm/ctn

Storage and Transportation Condition
-20°C~60°C(-4°F~140°F)
≤ 90% RH; 700~1060 hPa