



Automatic Wrist Blood Pressure Monitor

Model: S150



- Smart sense technology
- Movement Detection
- Hypertension Risk Indication
- Irregular Heartbeat Detection (IHB)
- Average of last 3 readings
- One-touch automatic operation
- 90 memories
- Carrying case



Smart Sense Technology



Movement Detection



Hypertension Risk Indication



Irregular Heartbeat Detection (IHB)



Average of last 3 readings



One-touch automatic operation



90 Memories



AAMI

This blood pressure monitor has passed the criteria of the European Society of Hypertension by achieving all requirements. This monitor further meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)



Cuff fits wrist circumference: 13.5~22cm (5.3"~8.7")

Smart Sense Technology

Smart Sense Technology utilizes an advanced air control system. The electronic deflation valve precisely controls the deflation rate and allows the monitor to detect and analyze pulse signals, adjusting cuff pressure during deflation. This significantly ensures more precise analysis, while reaching greater accuracy.

Model	Qty per carton	Carton volume
S150	36 pcs	0.047 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
21420 pcs	42876 pcs	49752 pcs

Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH, 700~1060 hPa