



HEALTH STATION	STL 501F (Deluxe)	STM 300C (Desktop)	STM 300F (Desktop)
	<b>Hardware Specifications</b>		
Size (incl. 4in1 rack)	(L) 720 x (D) 370 x (H) 894 mm	(L) 610 x (D) 315 x (H) 472 mm	
OS	Android 8.1		
CPU	1.5GHz		
RAM	2GB		
HHD	eMMC 16GB		
Network	LAN (Intel) x 1 - WLAN RTL8723BS 802.11 b/g/n 1T1R		
Bluetooth	4.0		
Speaker	2W		
Touch Screen	27" FHD Panel	15.6" FHD Panel	15.6" FHD Panel
Printer	●	●	●
4in1RACK	●	●	●
Camera	●	●	●
NFC Reader	●		●
Smart Card Reader	●		●
QR Code / Bar Code Scanner	●		●
	<b>Health Devices</b>		
X5 BT - Blood Pressure Monitor	●	●	●
HS200 - Blood Glucose Meter	●	●	●
SB210 - Pulse Oximeter	●	●	●
HC700 - Thermometer	●	●	●
WF262 - Body Fat Monitor	●	●	●
	<b>Software Specifications - Stand alone</b>		
Slideshow	●	●	●
Health Check Up	Blood Pressure, Blood Glucose, SpO <sub>2</sub> , Temperature, Weight, Body & Visceral Fat, BMI, BMR, Body Muscle		
Print Receipt	Check up result, Coupon		
Settings	Load slides, Set coupon content, Change banner		



MEDICAL KIT	STP 100C
	<b>Licensed Android System Software included</b>
<b>Recommended Tablet</b>	Samsung Galaxy Tab A7 Lite (T225) (not included)
	<b>Recommended Hardware Specifications for Tablets</b>
<b>OS</b>	Android 11.0 above
<b>CPU</b>	2GHz
<b>RAM</b>	2GB
<b>HHD</b>	32GB
<b>Network</b>	WiFi 802.11 b/g/n
<b>Bluetooth</b>	5.0
<b>Touch Screen</b>	<b>8.7" Panel (1340 x 800   1920 x 1080   1920 x 1200)</b>
<b>Camera</b>	For reading Barcodes
	<b>Health Devices included</b>
<b>X5 BT - Blood Pressure Monitor</b>	●
<b>HS200 - Blood Glucose Meter</b>	●
<b>SB210 - Pulse Oximeter</b>	●
<b>HC700 - Thermometer</b>	●
<b>WF262 - Body Fat Monitor</b>	●
	<b>Software Specifications - Stand alone</b>
<b>Health Check Up</b>	Blood Pressure, Blood Glucose, SpO <sub>2</sub> , Temperature, Weight, Body & Visceral Fat, BMI, BMR, Body Muscle
<b>Check Up Results</b>	By scanning the QR Code on the Tablet screen